




MEASUREMENTS

To refer on the prescription/Order form (page 1)


Rhinometer (red) , cephalometer (blue) , rule (black)






MEASUREMENTS on a frame

Refer to the prescription / order form page 2




MEASUREMENTS ON THE HEAD



SPHENOID WIDTH

Average « european » man :110-130mm /
woman :105-125mm




TEMPLE WIDTH

Av. M : 140-145/ Av. W :130-135





HEAD WIDTH



ROCHER WIDTH

NOSE ANGLES

FRONT ANGLE



Average in Europe : around 20°

SPLAY ANGLE



Average in Europe : between 20 and 25°